MARATHON PRODUCTS

c\temp-LCD™ Data Logger

Our Single-use Temperature Logger with Real-Time LCD Display for Instant Acceptance or Quarantine Decisions.

c\temp-LCD is the revolutionary, single-use, pre-programmed, electronic temperature data recorder specifically designed for measuring temperatures during the transportation of refrigerated vaccines, biological material, chemicals, pharmaceuticals, and for documenting the transportation of perishable commodities. The ctemp-LCD is NIST-Traceable and comes with a Certificate of Conformance.

Functions include:

- Internal Thermistor: accurately records time and temperatures for user selectable readings in Celsius or Fahrenheit.
- LCD display: Shows Current temperature, Max Temperature, Minimum Temperature, Run Time, Alarm Status, and Recording.

Small Form Factor: c\temp-lcd has an extended battery life and can endure rough handling. The internal software records the following: time mechanism, start time, sampling intervals, end description for special instructions, alarm values, recording status, measurement data, software version, serial and tracking number, model type and after the fact report generation options.

Visual Alarms: All our data loggers can be programmed with alarm limits. Should the temperature reach levels outside these limits, the LDC will display "ALARM." To start the recorder, simply press the start button!

USB Interface: To download recorded information, simply plug the USB reader station and cable into the standard USB port of your PC. MaxiThermal software, reader station and cable sold separately.

Cloud Computing

Temperature Logger Management

TRACEABILITY, HACCP COMPLIANCE


Marathon's Temperature data logger customers can share data and access their temperature and RH files remotely, anywhere, at anytime via the Internet.


Our cloud-based software provides compliance with HACCP rules and regulations for all of your perishable shipments.

CRYO STORE



MaxiThermal Software Highlights: Easy-to-use software provides clear graphs. Report Summary & Statistics File Manager, Note Attachment, Zoom, Elapsed Time or Date/Time View, Recorded View, Celsius or Fahrenheit. Data can be exported to other Windows programs such as spreadsheets, word-processing or databases. Available for Windows XP, Vista, 7, 8 and 10. An Enterprise Software version and 21CFR Part 11 Compliant Software are also available.





+31 (0) 513 41 59 65

for more information or to place an order.

www.cryostore.com

Cryo Store b.v. Handelswei 1, 8501 XJ Joure, The Netherlands VAT: NL819624032B01 • CoC: 01132538

c\temp-LCD Specifications SKU#7007

Memory: EEPROM, 8k Storage (approx. 8,000 data points)

Temperature Range: -35°C to $+50^{\circ}\text{C}$ (-31°F to +122°F)

Temperature Accuracy: $\pm 0.2^{\circ}\text{C}$ at 0° to 10°C , $\pm 1.0^{\circ}\text{C}$ at extremes

 $(\pm 0.36^{\circ}F$ at 32° to 50° F, $\pm 1.8^{\circ}F$ at extremes)

Measurement Intervals: Pre-programmed from 2 seconds to months.

Temperature Sensor: NIST-traceable Precision SMT Thermistor

Indicators: LCD Display. Start/Stop and Alarm

High and low alarm limits: ALARM set by temperature and time.

Computer Interface: USB Reader Station and Cable

Internal Firmware Provides: Internal clock, Start time, Interval, User description,

Alarm values, Recording status, Measurement data, Firmware version, Serial number, Model type

Power Source: One 3.0 Volt Lithium coin cell battery with a 2 year shelf life.

We recommend replacing the battery after 12 months of use.

Size: Length Width Depth

6.5 cm (2.56 in.) 3.4 cm (1.3 in.) 0.9 cm (0.35 in.)

Weight: 57.7 g (2.0 oz.)

Quantity: Minimum 20 units. Your order is shipped on the same

day if the order is placed before 2 p.m. PST.

Required Accessories: MaxiThermal Software

USB Reader Station and Cable

Serialization: This recorder is programmed at the factory with a unique

serial number for positive tracking.

Water-resistant Enclosure: The product package is hermetically sealed to an IP66 rating.

Recyclable: https://example.com/recorders are designed for single use, and we encourage recycling. Please call Marathon Products for complete disposal or recycling instructions.

